

# **HYBRID**

PORTABLE X-RAY SYSTEM

***Remember.  
Bring our care for you; Safety and Quality.***

**HYBRID S70**



**DIGIMED**

## Technical Specification

- Tube Voltage : 70 Kv
- Tube Current : 2 mA
- Focal Spot : 0.4 mm
- Time Range : 0.05 - 1 sec
- Total Filtration : Over 1.5 mmAl
- SSD : 200 mm
- Weight : 1.9 kg
- Dimension : 146 x 227 x 126 mm
- Power Source : 14.8V Rechargeable Battery

# DIGIMED



*Switchable  
Rectangular  
Emission cone*



**RELIABLE X-RAY SOURCE**  
Safety Requirements cleared



**EASY TO HANDLE**  
Best Grip on your palm  
Harmony with Patient



**ROTARY SWITCH**  
Easy and Fast Setting



**STRONGER SAFEGUARD**  
Internal Radiation Protection



**WIDE TFT COLOR LCD**  
More clearly visible

*Just right there  
you need it!*



*Distributed by,*

**ACTEON ANZ**

Suite 1, 13B Narabang Way Belrose NSW 2085

Office +61 2 9669 2292

**Toll Free 1300 Acteon** (1300 22 83 66)

Email [acteon.anz.admin@acteongroup.com](mailto:acteon.anz.admin@acteongroup.com)

*Dealers contact information:*