

# Purus 12 Litre B class Autoclave for Dental

Purus recognizes the importance to dentists of quality, reliability, and value for money.

Our range of autoclaves for dentists brings easy-to-use functionality, contemporary styling, technical quality, and the reassurance of warranty protection and service support.

This means that your practice can ensure the highest quality of sterilization standards which you and your patients deserve.

## 12 Litre B class for Dental

- European standard, with 3 times pre and drying vacuum
- Sterilizing Temp 121 °C 134 °C
- Residual humidity is <0.2% vacuum measurement reaches -0.70bar.
- Sterilizes:
  - miscellaneous packed or unpacked solid or hollow items
  - porous items
  - efficiently removes air from every hollow area of any instrument. Ensuring full sterilization.
- TGA approved ARTG 222444
- Conformance to:
  - prEN 10360/2/4
  - Australian standard AS/NZ 4815-2006
- CE compliant
- Total compliance for your sterilization needs



12 Litre B class Autoclave for Dental

## Product Code

purus\_AC-12B

## Product Description

### Chamber Size

- Diameter - 192mm
- Depth - 320mm

### External dimensions

- 56cm D x 37cm W x 36cm H

### Weight

- 43kg.

# 12 Litre B class Autoclave for Dental

## Installation

- Simply plug in and use – self-installed by customer
- Bowie Dick and Leak test Installed
- Accessories supplied:
  - printer
  - chamber rack
  - 3 trays
  - door spanner
  - tray tool
  - spare door seal
  - level

## Service

- Autoclave Calibration / Validation to be completed annually by a Purus Authorised Service Technician
- Quality and cost effective service supported by warranty
- Repair and Maintenance agreements available
- No hidden upkeep service requirements

## Warranty

- 2 year warranty on Autoclave
- 12 month warranty on printer
- Optional extended warranty available.

## Finance

We have a panel of finance providers (including Medfin, Westpac Health, and Investec) and will introduce you to them for a competitive financing package.

### Autoclave calibration graph



### Contact Details

Contact 1300 Purus 9  
Unit 4 / 56 Barclay Road, Derrimut VIC 3030  
**p:** +61 3 8358 5557  
**f:** +61 3 8358 5545  
**e:** [info@purushealth.com.au](mailto:info@purushealth.com.au)  
[www.purushealth.com.au](http://www.purushealth.com.au)

## Call Purus Now

- To discuss your Autoclave needs
- For best pricing