

VibraKleen E2 Ingredient List



The Natural Power of enzymes

VibraKleen E2 is an enzymatic ultrasonic cleaning solution that is available in easy-to-use tablet dispensing.

VibraKleen E2 tablets are convenient and effective for:

- Pre-soak of instruments
- Evacuation cleaner

VibraKleen E2 is available in boxes of 20 or Clinic Pack of 80.

Product Name: VibraKleen E2 Tablets

Product Code:

Ingredient	CAS No.	Function
Subtilisin	9014-01-1	Protease enzyme, buffer
Sodium Carbonate	497-19-8	Alkalizer, buffer
Citric Acid	77-92-9	Preservative, effervescing agent, pH control agent
Sodium Benzoate	532-32-1	Preservative
Polyethylene Glycol	25322-68-3	Filler
Spearmint Flavor Natural NV-9466	Not Available	Fragrance ingredient
FD&C Blue#1	3844-45-9	Colorant

Giving a hand to oral health.



KULZER
MITSUI CHEMICALS GROUP